

FORMFIT® WALKER AIR

The Formfit Walker Air provides immobilisation for patients suffering from soft tissue injury; grades 2 & 3 strains and sprains, and stable fractures. The patented pneumatic system is designed to increase stability while decreasing pain and swelling. The intuitive pump and clearly marked decompression valve make air/pressure adjustments easy.

INDICATIONS FOR USE

For foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Soft tissue injuries
- During postoperative or rehab phase

FEATURES AND BENEFITS

- Available in low top or standard height
- Integrated pneumatic system
- Low, wide rocker bottom promotes a stable, natural gait
- Shock absorbing sole for reduced impact at heel strike
- Breathable, seamless liner

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

Size	Shoe Size UK	Shoe Size EU	Standard	Low Top
X-Small	1.5 - 3.5	33 - 35.5	A-W0200BLK	A-W0300BLK
Small	4 - 6.5	36 - 39	A-W0400BLK	A-W0500BLK
Medium	7 - 10	39.5 - 44	A-W0600BLK	A-W0700BLK
Large	10.5 - 12	44.5 - 46.5	A-W0800BLK	A-W0900BLK
X-Large	12.5 - 14.5	47 - 50	A-W1000BLK	A-W1100BLK

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Standard	Low Top
X-Small	12349	12356
Small	10675	12357
Medium	10678	12358
Large	10681	12359
X-Large	12347	12360

OTHER ACCESSORIES

Description	Part#
10cm softgood width extension	WKSG-1BLK
Tibia and calf plates for use when additional splinting is required	WKSP100
Hot/Cold therapy insert for safe and reusable cryotherapy	WKGEL100
Heel wedge	AWK0005
Adhesive insole for heel wedge	WHW-W0107